

Manufacturer Disclosure Statement for Medical Device Security – MDS²

DEVICE DESCRIPTION

Device Category Ultrasound	Manufacturer Butterfly Network, Inc.	Document ID 950-20007-00 v05	Document Release Date November, 2019
Device Model iQ	Software Revision 1.13.0	Software Release Date October, 2019	
Manufacturer or Representative Contact Information	Company Name Butterfly Network, Inc. Representative Name/Position Allan Bottemiller Medical Informatics Manager	Manufacturer Contact Information Phone: (203) 458-7100 Email: support@butterflynetwork.com Web: https://butterflynetwork.com/contact	

Intended use of device in network-connected environment:
 Butterfly Network, Inc. (BNI) develops and manufactures point of care ultrasound (POCUS) products. Butterfly iQ, the POCUS device, has been cleared by FDA as a Class II medical device. Ultrasound studies are conducted using the iQ device and scans are visualized and controlled from the Butterfly mobile app. The Butterfly Cloud allows for the storage and transmission of ultrasound images, integration with hospital DICOM end points (PACS/VNA), EMR via HL7, and facilitates cross-team HIPAA-compliant communication. The Butterfly Cloud leverages best in class SaaS technologies to ensure maximum uptime and security (e.g. Amazon Web Services ‘AWS’, Aptible, Auth0). In transit, all data is protected with TLS 1.2 encryption standard. At rest, all data is protected by AES 256-bit encryption. All products are developed under the Butterfly quality management system (QMS). Butterfly Network is HIPAA and HITECH compliant with SOC II certification.

MANAGEMENT OF PRIVATE DATA

Refer to Section 2.3.2 of this standard for the proper interpretation of information requested in this form.			Yes, No, N/A, or See Note	Note #
A	Can this device display, transmit, or maintain private data (including electronic Protected Health Information [ePHI])?		Yes	1
B	Types of private data elements that can be maintained by the device :			
	B.1 Demographic (e.g., name, address, location, unique identification number)?		Yes	
	B.2 Medical record (e.g., medical record #, account #, test or treatment date, device identification number)?		Yes	
	B.3 Diagnostic/therapeutic (e.g., photo/radiograph, test results, or physiologic data with identifying characteristics)?		Yes	
	B.4 Open, unstructured text entered by device user/operator ?		Yes	
	B.5 Biometric data ?		No	
	B.6 Personal financial information?		No	
C	Maintaining private data - Can the device :			
	C.1 Maintain private data temporarily in volatile memory (i.e., until cleared by power-off or reset)?		Yes	2
	C.2 Store private data persistently on local media?		See Note	3
	C.3 Import/export private data with other systems?		Yes	4
	C.4 Maintain private data during power service interruptions?		Yes	
D	Mechanisms used for the transmitting, importing/exporting of private data – Can the device :			
	D.1 Display private data (e.g., video display, etc.)?		Yes	
	D.2 Generate hardcopy reports or images containing private data ?		No	
	D.3 Retrieve private data from or record private data to removable media (e.g., disk, DVD, CD-ROM, tape, CF/SD card, memory stick, etc.)?		No	
	D.4 Transmit/receive or import/export private data via dedicated cable connection (e.g., IEEE 1073, serial port, USB, FireWire, etc.)?		No	
	D.5 Transmit/receive private data via a wired network connection (e.g., LAN, WAN, VPN, intranet, Internet, etc.)?		No	
	D.6 Transmit/receive private data via an integrated wireless network connection (e.g., WiFi, Bluetooth, infrared, etc.)?		Yes	
	D.7 Import private data via scanning?		Yes	
	D.8 Other?		-----	---

Management of Private Data notes:

1. The Butterfly iQ provides the ability for users to manually enter PHI or pull data from a DICOM Modality Worklist.
2. Ultrasound exams that may include PHI are cached within the Butterfly mobile app until successfully uploaded to the Butterfly Cloud.
3. Users are required to access the internet on the device at least once every 30 days to validate the Butterfly mobile app, iQ device and account registration. If there is no access to the internet from the mobile device for more than 30 days, a notice will persist and no other buttons will be accessible. At this point, the user is blocked from accessing any data within the Butterfly mobile app until they go online.
4. The Butterfly Cloud utilizes HTTPS with TLS 1.2 encryption for hospital system connections, such as PACS for archiving and Modality Worklist to pull patient demographics to the Butterfly iQ mobile app.

Device Category Ultrasound	Manufacturer Butterfly Network, Inc.	Document ID 950-20007-00 v05	Document Release Date November, 2019
Device Model iQ	Software Revision 1.13.0	Software Release Date October, 2019	

SECURITY CAPABILITIES

Refer to Section 2.3.2 of this standard for the proper interpretation of information requested in this form. Yes, No,
N/A, or
See Note Note #

1 AUTOMATIC LOGOFF (ALOF)

The **device's** ability to prevent access and misuse by unauthorized **users** if **device** is left idle for a period of time.

1-1	Can the device be configured to force reauthorization of logged-in user(s) after a predetermined length of inactivity (e.g., auto-logoff, session lock, password protected screen saver)?	Yes	1
-----	---	-----	---

1-1.1	Is the length of inactivity time before auto-logoff/screen lock user or administrator configurable? (Indicate time [fixed or configurable range] in notes.)	Yes	2
-------	--	-----	---

1-1.2	Can auto-logoff/screen lock be manually invoked (e.g., via a shortcut key or proximity sensor, etc.) by the user ?	Yes	1, 2
-------	---	-----	------

ALOF notes:

1. The inactivity timer is set within the mobile device configuration. The Butterfly mobile app forces users to apply a security PIN, password, etc. to their mobile device. This can be enforced by the use of a Mobile Device Management (MDM) software as well. Removing the device lock will result in the user being automatically logged out of the Butterfly app, which will permanently delete any iQ data that has not been uploaded to the Butterfly Cloud. The user will be unable to log back into the Butterfly mobile app until a device lock has been re-enabled.
2. The Butterfly Cloud provides a customizable inactivity timer that can be set by an administrator between 15 minutes and 24 hours.

2 AUDIT CONTROLS (AUDT)

The ability to reliably audit activity on the **device**.

2-1	Can the medical device create an audit trail ?	Yes	1
-----	--	-----	---

2-2	Indicate which of the following events are recorded in the audit log:		
-----	---	--	--

2-2.1	Login/logout	Yes	___
-------	--------------	-----	-----

2-2.2	Display/presentation of data	See Note	2
-------	------------------------------	----------	---

2-2.3	Creation/modification/deletion of data	Yes	___
-------	--	-----	-----

2-2.4	Import/export of data from removable media	N/A	___
-------	---	-----	-----

2-2.5	Receipt/transmission of data from/to external (e.g., network) connection	Yes	___
-------	--	-----	-----

2-2.5.1	Remote service activity	N/A	3
---------	--------------------------------	-----	---

2-2.6	Other events? (describe in the notes section)	See Note	4
-------	---	----------	---

2-3	Indicate what information is used to identify individual events recorded in the audit log:		
-----	--	--	--

2-3.1	User ID	Yes	___
-------	----------------	-----	-----

2-3.2	Date/time	Yes	___
-------	-----------	-----	-----

AUDT notes:

1. Any data accessed on the Butterfly Cloud and transmitted to the Butterfly Cloud is logged.
2. User activity is logged, including which ultrasounds studies have been accessed.
3. There is no remote access to the device through the Butterfly mobile app.
4. All user activity (logging in, logging out, screens viewed, date/time of new uploads, deletions, etc. on studies, captures and archives) is logged.

3 AUTHORIZATION (AUTH)

The ability of the device to determine the authorization of users.

3-1	Can the device prevent access to unauthorized users through user login requirements or other mechanism?	Yes	1, 2
-----	--	-----	------

3-2	Can users be assigned different privilege levels within an application based on 'roles' (e.g., guests, regular users , power users , administrators, etc.)?	Yes	___
-----	--	-----	-----

3-3	Can the device owner/operator obtain unrestricted administrative privileges (e.g., access operating system or application via local root or admin account)?	No	___
-----	--	----	-----

AUTH notes:

1. The Butterfly mobile app requires a device lock (password, pin, etc.) to be enabled. SSO/SAML integration can be implemented to enforce user login and password requirements.
2. Mobile device restrictions can be enforced through the use of the Enterprise Mobility Management (EMM) feature on the Butterfly Cloud to prevent users from logging into the Butterfly mobile app on personal devices.

Device Category	Manufacturer	Document ID	Document Release Date
Ultrasound	Butterfly Network, Inc.	950-20007-00 v05	November, 2019
Device Model	Software Revision	Software Release Date	
iQ	1.13.0	October, 2019	
Refer to Section 2.3.2 of this standard for the proper interpretation of information requested in this form.			Yes, No, N/A, or See Note
Note #			
4	CONFIGURATION OF SECURITY FEATURES (CNFS)		
	The ability to configure/re-configure device security capabilities to meet users' needs.		
4-1	Can the device owner/operator reconfigure product security capabilities ?		Yes 1
CNFS notes:	1. The device owner/operator can only configure security capabilities local to the mobile device, such as Auto Lock and password. User access to the Butterfly mobile app and data stored on the Butterfly Cloud is controlled by a designated Butterfly Cloud administrator.		
5	CYBER SECURITY PRODUCT UPGRADES (CSUP)		
	The ability of on-site service staff, remote service staff, or authorized customer staff to install/upgrade device's security patches.		
5-1	Can relevant OS and device security patches be applied to the device as they become available?		Yes ___
	5-1.1	Can security patches or other software be installed remotely?	See Note 1
CSUP notes:	1. The Butterfly mobile app does not impact mobile device security patches from being applied/installed. Butterfly Network security patches are delivered through the Butterfly mobile app.		
6	HEALTH DATA DE-IDENTIFICATION (DIDT)		
	The ability of the device to directly remove information that allows identification of a person.		
6-1	Does the device provide an integral capability to de-identify private data ?		Yes 1
DIDT notes:	1. Users have the ability to share de-identified study links from the Butterfly Cloud and Butterfly mobile app.		
7	DATA BACKUP AND DISASTER RECOVERY (DTBK)		
	The ability to recover after damage or destruction of device data, hardware, or software.		
7-1	Does the device have an integral data backup capability (i.e., backup to remote storage or removable media such as tape, disk)?		Yes 1
DTBK notes:	1. Ultrasound exams acquired on the Butterfly iQ are securely uploaded to the Butterfly Cloud, which is backed up daily and prior to updates.		
8	EMERGENCY ACCESS (EMRG)		
	The ability of device users to access private data in case of an emergency situation that requires immediate access to stored private data .		
8-1	Does the device incorporate an emergency access ("break-glass") feature?		No ___
EMRG notes:	N/A		
9	HEALTH DATA INTEGRITY AND AUTHENTICITY (IGAU)		
	How the device ensures that data processed by the device has not been altered or destroyed in an unauthorized manner and is from the originator.		
9-1	Does the device ensure the integrity of stored data with implicit or explicit error detection/correction technology?		Yes 1
IGAU notes:	1. Data at rest is encrypted with AES 256-bit and data in transit is encrypted with TLS 1.2.		

Device Category Ultrasound	Manufacturer Butterfly Network, Inc.	Document ID 950-20007-00 v05	Document Release Date November, 2019
Device Model iQ	Software Revision 1.13.0	Software Release Date October, 2019	

Refer to Section 2.3.2 of this standard for the proper interpretation of information requested in this form.

Yes, No,
N/A, or
See Note

Note #

10 MALWARE DETECTION/PROTECTION (MLDP)			
The ability of the device to effectively prevent, detect and remove malicious software (malware).			
10-1	Does the device support the use of anti-malware software (or other anti-malware mechanism)?		See Note 1
	10-1.1	Can the user independently re-configure anti-malware settings?	Yes —
	10-1.2	Does notification of malware detection occur in the device user interface?	See Note 1
	10-1.3	Can only manufacturer-authorized persons repair systems when malware has been detected?	N/A 1
10-2	Can the device owner install or update anti-virus software ?		Yes —
10-3	Can the device owner/ operator (technically/physically) update virus definitions on manufacturer-installed anti-virus software ?		N/A —
MLDP notes: 1. Installing anti-malware software on the mobile device will not impact the Butterfly mobile app. Anti-malware notifications may appear on the device user interface, but not from within the Butterfly mobile app.			

11 NODE AUTHENTICATION (NAUT)			
The ability of the device to authenticate communication partners/nodes.			
11-1	Does the device provide/support any means of node authentication that assures both the sender and the recipient of data are known to each other and are authorized to receive transferred information?		Yes 1
NAUT notes: 1. Transmission of data between the Butterfly mobile app and Butterfly Cloud solution is performed using HTTPS with TLS 1.2 encryption. DICOM Store and Modality Worklist (MWL) connections are performed via Butterfly Link software using HTTPS with TLS 1.2 encryption or DICOM-TLS (TLS 1.2).			

12 PERSON AUTHENTICATION (PAUT)			
Ability of the device to authenticate users			
12-1	Does the device support user/operator -specific username(s) and password(s) for at least one user ?		Yes 1
	12-1.1	Does the device support unique user/operator -specific IDs and passwords for multiple users?	Yes —
12-2	Can the device be configured to authenticate users through an external authentication service (e.g., MS Active Directory, NDS, LDAP, etc.)?		Yes 2
12-3	Can the device be configured to lock out a user after a certain number of unsuccessful logon attempts?		Yes —
12-4	Can default passwords be changed at/prior to installation?		N/A 3
12-5	Are any shared user IDs used in this system?		No —
12-6	Can the device be configured to enforce creation of user account passwords that meet established complexity rules?		Yes 2
12-7	Can the device be configured so that account passwords expire periodically?		Yes 2
PAUT notes: 1. Butterfly Network user accounts are managed by a designated Butterfly Cloud administrator. 2. SSO/SAML integration is available to enforce corporate user login and password requirements. 3. There are no default passwords.			

13 PHYSICAL LOCKS (PLOK)			
Physical locks can prevent unauthorized users with physical access to the device from compromising the integrity and confidentiality of private data stored on the device or on removable media .			
13-1	Are all device components maintaining private data (other than removable media) physically secure (i.e., cannot remove without tools)?		Yes —
PLOK notes: N/A			

Device Category	Manufacturer	Document ID	Document Release Date	
Ultrasound	Butterfly Network, Inc.	950-20007-00 v05	November, 2019	
Device Model	Software Revision	Software Release Date		
iQ	1.13.0	October, 2019		
Refer to Section 2.3.2 of this standard for the proper interpretation of information requested in this form.			Yes, No, N/A, or See Note	Note #
14	ROADMAP FOR THIRD PARTY COMPONENTS IN DEVICE LIFE CYCLE (RDMP)			
	Manufacturer's plans for security support of 3rd party components within device life cycle.			
14-1	In the notes section, list the provided or required (separately purchased and/or delivered) operating system(s) - including version number(s).	See Note	1, 2	
14-2	Is a list of other third party applications provided by the manufacturer available?	Yes	___	
RDMP notes:	<p>1. A list of supported mobile devices can be found at https://www.butterflynetwork.com/specs</p> <p>2. The Butterfly Cloud solution is compatible with the latest two versions of Google Chrome, Mozilla Firefox, Safari and Microsoft Edge.</p>			
15	SYSTEM AND APPLICATION HARDENING (SAHD)			
	The device's resistance to cyber attacks and malware .			
15-1	Does the device employ any hardening measures? Please indicate in the notes the level of conformance to any industry-recognized hardening standards.	See Note	1	
15-2	Does the device employ any mechanism (e.g., release-specific hash key, checksums, etc.) to ensure the installed program/update is the manufacturer-authorized program or software update?	Yes	2	
15-3	Does the device have external communication capability (e.g., network, modem, etc.)?	Yes	3	
15-4	Does the file system allow the implementation of file-level access controls (e.g., New Technology File System (NTFS) for MS Windows platforms)?	No	___	
15-5	Are all accounts which are not required for the intended use of the device disabled or deleted, for both users and applications?	Yes	4	
15-6	Are all shared resources (e.g., file shares) which are not required for the intended use of the device , disabled?	No	5	
15-7	Are all communication ports which are not required for the intended use of the device closed/disabled?	Yes	___	
15-8	Are all services (e.g., telnet, file transfer protocol [FTP], internet information server [IIS], etc.), which are not required for the intended use of the device deleted/disabled?	Yes	___	
15-9	Are all applications (COTS applications as well as OS-included applications, e.g., MS Internet Explorer, etc.) which are not required for the intended use of the device deleted/disabled?	No	5	
15-10	Can the device boot from uncontrolled or removable media (i.e., a source other than an internal drive or memory component)?	No	___	
15-11	Can software or hardware not authorized by the device manufacturer be installed on the device without the use of tools?	No	___	
SAHD notes:	<p>1. Butterfly Network follows ISO 27001 standards and is built under the Butterfly QMS.</p> <p>2. Software and firmware updates are only available through the Butterfly mobile app.</p> <p>3. Cellular and WiFi are supported.</p> <p>4. There are no default accounts. Access to the Butterfly app and Butterfly Cloud is restricted to active user accounts, which are managed by a designated Butterfly Cloud administrator for each organization.</p> <p>5. The Butterfly mobile app allows the mobile device to function normally, however, the app itself only utilizes specified ports and communication protocols when in use.</p>			
16	SECURITY GUIDANCE (SGUD)			
	The availability of security guidance for operator and administrator of the system and manufacturer sales and service.			
16-1	Are security-related features documented for the device user ?	Yes	___	
16-2	Are instructions available for device/media sanitization (i.e., instructions for how to achieve the permanent deletion of personal or other sensitive data)?	Yes	1	
SGUD notes:	<p>1. Data is only cached securely within the Butterfly mobile app until uploaded to the Butterfly Cloud. The Butterfly mobile app requires the device to have a lock mechanism enabled before the app can be used. Removing the device lock will result in the user being automatically logged out of the Butterfly app, which will permanently delete any iQ data that has not been uploaded to the Butterfly Cloud. The user will be unable to log back into the Butterfly mobile app until a device lock has been re-enabled.</p>			

Device Category	Manufacturer	Document ID	Document Release Date	
Ultrasound	Butterfly Network, Inc.	950-20007-00 v05	November, 2019	
Device Model	Software Revision	Software Release Date		
iQ	1.13.0	October, 2019		
Refer to Section 2.3.2 of this standard for the proper interpretation of information requested in this form.			Yes, No, N/A, or See Note	Note #
17	HEALTH DATA STORAGE CONFIDENTIALITY (STCF)			
	The ability of the device to ensure unauthorized access does not compromise the integrity and confidentiality of private data stored on device or removable media .			
17-1	Can the device encrypt data at rest?		Yes	1
STCF notes:	1. Data cached within the Butterfly mobile app is encrypted with AES 256-bit.			
18	TRANSMISSION CONFIDENTIALITY (TXCF)			
	The ability of the device to ensure the confidentiality of transmitted private data .			
18-1	Can private data be transmitted only via a point-to-point dedicated cable?		No	—
18-2	Is private data encrypted prior to transmission via a network or removable media ? (If yes, indicate in the notes which encryption standard is implemented.)		Yes	1, 2
18-3	Is private data transmission restricted to a fixed list of network destinations?		Yes	2
TXCF notes:	1. Transmission of data between Butterfly mobile app and the Butterfly Cloud solution is performed using HTTPS with TLS 1.2 encryption. 2. DICOM Store and Modality Worklist (MWL) connections are performed via Butterfly Link software using HTTPS with TLS 1.2 encryption or DICOM-TLS (TLS 1.2).			
19	TRANSMISSION INTEGRITY (TXIG)			
	The ability of the device to ensure the integrity of transmitted private data .			
19-1	Does the device support any mechanism intended to ensure data is not modified during transmission? (If yes, describe in the notes section how this is achieved.)		Yes	1
TXIG notes:	1. Transmission of data between Butterfly mobile app and the Butterfly Cloud solution is performed using HTTPS with TLS 1.2 encryption.			
20	OTHER SECURITY CONSIDERATIONS (OTHR)			
	Additional security considerations/notes regarding medical device security.			
20-1	Can the device be serviced remotely?		No	—
20-2	Can the device restrict remote access to/from specified devices or users or network locations (e.g., specific IP addresses)?		N/A	—
	20-2.1 Can the device be configured to require the local user to accept or initiate remote access?		N/A	—
OTHR notes:	N/A			